

WOUND WASH

- 100 ML -

WOUND WASH

XX Wound Wash provides a gentle spray for irrigation and cleansing of wounds. It cleans the wound and, in that way, supports wound healing without the need of touching the wound.

Intended use
Sterile isotonic saline solution for rinsing/cleansing of wounds.

Target user
For all ages. To be used for self-care, at institutions or workplaces without technical experience or training.

Directions
Spray a small quantity through the nozzle before each use.
Spray the area to be treated from a distance of 10-20 cm, working across the wound to thoroughly cleanse and irrigate the area.

Composition
XXX

Storage
Store between 15°C (59°F) and 30°C (86°F) for comfortable use.
Storage/Transport temperature: -20°C (-4°F) to 40°C (104°F).

Disposal
Follow local regulations for safe disposal of pressurized containers.

Precautions

- Avoid direct contact between the nozzle and the surface of the wound.
- Seek medical attention if condition persists or gets worse.
- If any serious incident occurs in relation to the device, report to the local competent authority and the manufacturer

Warning
Pressurised container: May burst if heated.

- Do not pierce or burn, even after use.
- Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
- Do not expose to temperatures exceeding 50°C/122°F.
- Protect from sunlight.
- Keep out of reach of children.

3 STERILE R REF 2001-X

Aurena Laboratories AB
Fjärriksvägen 22
SE-653 50 Karlstad, Sweden

CE MD 175
2862

(01)XXXXXXXXXXXXC

LOT (10) See bottom of can
(17)

Distributor: XXXX

V4, 2022-12-01AP