

XX

Brim. volume (brimful capacity, total capacity)



Caution.



Distributor symbol. The symbol is accompanied by the distributor's name and address.

3

Inverted Epsilon.



Sterilized using irradiation



Reference number. This symbol must be followed by the reference number for the Aurena product family.



Manufacturer, CE mark, name and full address of manufacturer.  
The CE mark shall be placed next to the manufacturer symbol with the address. CE mark must have a minimum height of 5 mm.



2862



Medical Device.

(01) XXXXXXXXXXXXXC

(01) Article code, shall be placed above the 'LOT' and 'EXP' symbols.  
The UDI symbol shall be placed in front of (01) if multiple data carriers are present on the label.



The batch code symbol 'LOT' (10) and the use-by date symbol 'EXP' (17) shall be included in the marking together with a reference to; 'See bottom of can'. The same applies for the carton box; 'See bottom of box' and in the leaflet; 'See bottom of can/box'. E.g.: