

Wound Wash

Sterile saline solution for cleansing/irrigation of wounds that mechanically removes foreign bodies and reduces the number of pathogens in the wound



FACT BOX

Wound Wash

Market segments:

Healthcare institutions, community care, OTC

Composition: Sterile 0.9% NaCl solution or Sterile Ringer solution.

Product approvals: Registered according to MDD (Medical Device Directive, 93/42/EEC, Class IIb)

Shelf life: 3 years from production date.

Packaging: Standard can sizes are: 100, 125, 150, 200 and 250ml. Individual carton boxes, with or without leaflets are also available.

Presentations: Both sticker labels and preprinted cans can be used. Labels can be provided with a tamper seal. Preprinted cans are recommended for quantities of 25,000+ per order.

Wound healing is supported by this saline spray which mechanically removes foreign bodies, reduces pathogens, as well as blood, clot and dead tissue and moistening the tissue in the process. All this helps to reduce infection in the wound, which is one of the most frequent complications of wound healing.

Since the propellant, compressed air or nitrogen, is not flammable, the spray is classed as a non-flammable aerosol, which offers advantages regarding safety, transportation and storage.

Description:

Sterile isotonic saline solution for rinsing/cleansing of wounds. For all ages.

Claims:

- The spray may be used on acute and chronic wounds and pin sites to be cleansed
- The spray fluid softens debris and eschar to enable removal
- The spray provides a gentle way of removing debris, eschar, pus and slough from a wound cavity
- Support wound healing
- The spray works in any angle
- Can be used until last drop
- Hypoallergenic
- Sterile until last drop



AURENA
LABORATORIES

Fjärrviksvägen 22
653 50 Karlstad
Sweden

WWW.AURENALABS.COM

Version 2020-12-17