

# Ear Spray

Ear cleansing spray that prevents and removes the buildup of earwax



## FACT BOX

Ear Spray

**Market segments:** OTC

**Composition:** Sterile isotonic seawater solution or alternatively: sterile isotonic saline solution

**Product approvals:**  
Registered according to MDR  
(Medical Device Regulation  
2017/745, Class IIa)

**Shelf life:** 3 years from  
production date.

**Packaging:** Standard can sizes  
are: 50, 100, 125, 150, 200  
and 250 ml. Individual carton  
boxes, with or without leaflets  
are also available.  
Pre-printed can.

The ear cleansing spray from Aurena Laboratories is a medical device that cleans the ear and removes the buildup of earwax in a safe and convenient way.

EarSpray help cleanse the ear from wax or debris. The spray also clean ears without harming the ear together with preventing the buildup of more earwax and thus helps to maintain clean and healthy ears. The spray is easy to use and can be sprayed in any angle until the last drop. The formula contains no harsh chemicals and is preservative free.

Since the propellant, compressed air or nitrogen, is not flammable, the spray is classed as a non-flammable aerosol, which offers advantages regarding safety, transportation and storage.

### Intended use:

For preventive control of wax accumulation in the external ear canal.  
Can be used by adults and children from 30 months of age.

### Clinical claims:

- Prevention of build-up of ear wax
- Removal of ear wax
- Rinsing away ear wax, without pressure to the ear drum
- Helps maintain healthy ears

### Product characteristics :

- Sterilized
- Can be used at any angle
- Sterile to the last drop
- Preservative-free



**AURENA**  
LABORATORIES

Fjärrviksvägen 22  
653 50 Karlstad  
Sweden

[WWW.AURENALABS.COM](http://WWW.AURENALABS.COM)

V3. 2023-06-08AP