

Nasal Spray

Cleanse the nasal cavity and restore moisture to the nasal passages



FACT BOX

Nasal Sprays

Market segments: OTC

Composition: See separate table

Approved age: See separate table

Product approvals: Registered according to MDR (Medical Device Regulation 2017/745, Class IIa)

Shelf life: 3 years from production date.

Packaging: Standard can sizes are: 50, 100, 125, 150 and 200ml. Pre-printed cans. If a carton box, with or without leaflet, is required, they are available for can sizes up to 150 ml.

Aurena Laboratories have developed a broad range of nasal spray products based on sterile saline solution, sterile sea salt solution and natural sea water. The nasal sprays help to restore the moisture to nasal passages irritated by colds or allergies, dryness, pollution and overuse of nasal decongestants.

The seawater/saline solutions is a non-habit forming preservative-free nasal irrigation spray. The spray keeps the nasal passages clean, moisturized, helps relieve nasal congestion and clears the nasal cavity prior to administration of medicinal products. Can be used as frequent as needed.

Intended use:

To clear and moisturize the nasal mucosa/cavity.

Clinical claims:

- Relieves nasal congestion
- Alleviates allergy-related symptoms such as nasal congestion
- Cleanse the nasal cavity from allergens, irritants and bacteria
- Enables removal of obstructing matter
- Cleanse the nasal cavity prior to applying medical products
- Suitable for daily use
- Moisturizes the nasal cavity
- Reduces biofilms

Product characteristics:

- Sterile to the last drop
- Preservative-free
- Drug free
- Can be used at any angle
- Seawater contains minerals and trace elements
- Addition of aloe vera and/or chamomile to give a more pleasant cleaning in terms of feeling, scent and taste



Fjärrviksvägen 22
653 50 Karlstad
Sweden

